


### Pre-fabricated Wrist Extension Brace

*How to measure for a Pre-fabricated Wrist Brace*

- See diagram to check proper positioning for measurements.

**Level A** - Circumference at forearm band/strap location.

**Level B** - M-L of hand across MP joints.



**SMALL**

	Circumference	Diameter (M-L)
LEVEL A	6 1/2" - 9"	
LEVEL B		2 1/2" - 3 1/4"

**MEDIUM**

LEVEL A	9" - 11"	
LEVEL B		3 1/4" - 4 1/4"